

## MYKITA DIY PROTECTIVE SHIELD: ADVANCED


Designed as a practical contribution to current health initiatives, the below instructions are intended to provide a temporary protective solution in the absence of certified medical supplies (\* please refer to our full disclaimer below).




1  
You will need: kitchen sponge, elastic band, scissors, stapler, double-sided tape and an A4 clear plastic sheet



2  
Cut the sponge in half lengthways.



3  
Apply a strip of double-sided sticky tape to the long side of the plastic sheet as shown.



4  
Mount both sponge pieces to the tape as shown.



5  
Cut a piece of elastic to suit your head measurement. Fold the top corners of the sheet inwards, place the elastic in-between and secure with staples.



6  
Fold the bottom part of the sheet inwards and staple the pointed edge together.



7  
Position the soft part of the sponge on the forehead and secure via the elastic.



8  
Your DIY protective face shield is ready.

\*Disclaimer of liability: The manual for the DIY protective face shields is intended to provide temporary solutions when medical protective equipment is not available. The shields are designed to limit exposure of liquids touching sensitive areas in the face; they offer no protection against airborne viruses. They are not intended for children or persons in need of supervision and must be kept out of their hands. Clean and virus-free materials and tools must be used; the official health & safety guidelines by the World Health Organisation must be observed at all times. The DIY concepts are neither tested nor certified. As such, MYKITA explicitly rejects any liability for the effectiveness, manufacture or proper use of the protective shields. Any claim for damages against MYKITA due to injury to life, body or health – irrespective of the user or respective counterparty – is prohibited. The manufacture and use of the DIY protective shields is the sole responsibility of the user.