

MYKITA DIY PROTECTIVE FACE SHIELD: CAP

Designed as a practical contribution to current health initiatives, the below instructions are intended to provide a temporary protective solution in the absence of certified medical supplies (* please refer to our full disclaimer below).




1 You will need: a cap, duct tape, and an A4 clear plastic sheet.



2 Cut a piece of duct tape large enough to cover the gap at the back of the cap.



3 Secure the tape on the inside of the cap to cover the gap.



4 Mount the long side of the plastic sheet with duct tape, leaving half the sticky side exposed.



5 Mount the plastic sheet along the back of the cap using the sticky edge of the duct tape.



6 Press the sticky sides of the duct tape in the gap together for a firm hold.



7 Wear with clear plastic sheet in front to limit exposure to viral droplets.



8 Your DIY Protective Face Shield is ready.

*Disclaimer of liability: The manual for the DIY protective face shields is intended to provide temporary solutions when medical protective equipment is not available. The shields are designed to limit exposure of liquids touching sensitive areas in the face; they offer no protection against airborne viruses. They are not intended for children or persons in need of supervision and must be kept out of their hands. Clean and virus-free materials and tools must be used; the official health & safety guidelines by the World Health Organisation must be observed at all times. The DIY concepts are neither tested nor certified. As such, MYKITA explicitly rejects any liability for the effectiveness, manufacture or proper use of the protective shields. Any claim for damages against MYKITA due to injury to life, body or health – irrespective of the user or respective counterparty – is prohibited. The manufacture and use of the DIY protective shields is the sole responsibility of the user.