

## MYKITA DIY PROTECTIVE FACE SHIELD: TOWEL

Designed as a practical contribution to current health initiatives, the below instructions are intended to provide a temporary protective solution in the absence of certified medical supplies (\* please refer to our full disclaimer below).




1  
You will need: a dish towel or rectangular piece of fabric and an A4 clear plastic sheet.



2  
Lay the towel in a horizontal position and fold the top and bottom corners as shown.




3  
You want to create an elongated rectangle.



4  
Again, fold the long edges from top and bottom to meet in the middle as shown.



5  
Fold the towel over the plastic sheet as shown to secure.



6  
Make a loose knot of the towel ends.



7  
Position the shield on the forehead and tighten the knot to secure.



8  
Your DIY protective face shield is ready.

\*Disclaimer of liability: The manual for the DIY protective face shields is intended to provide temporary solutions when medical protective equipment is not available. The shields are designed to limit exposure of liquids touching sensitive areas in the face; they offer no protection against airborne viruses. They are not intended for children or persons in need of supervision and must be kept out of their hands. Clean and virus-free materials and tools must be used; the official health & safety guidelines by the World Health Organisation must be observed at all times. The DIY concepts are neither tested nor certified. As such, MYKITA explicitly rejects any liability for the effectiveness, manufacture or proper use of the protective shields. Any claim for damages against MYKITA due to injury to life, body or health – irrespective of the user or respective counterparty – is prohibited. The manufacture and use of the DIY protective shields is the sole responsibility of the user.